
E940 med.

112-150

How to buy

- **1. Measure** your client to determine the size you need. (See dimension chart for help.)
- **2. Select** the appropriate size and color.
- **3. Choose** the accessories you need.



Key user dimension: height

Select the appropriate Prone Stander by the user's overall height. Choose the model that allows for growth.

E930 small

Important: User's weight must not exceed the maximum working load.

Item dimensions (cm)	E930 small	E940 med.	E950 large
Board length	60 (41 without kneeboard)	76	99
Board width (top/bottom)	32/32	28/41	36/48
Top of board to footboard	51-89	84-110	112 – 145
Height when horizontal	75	75	80
Height when vertical	98	116	155
Distance between laterals	20-41	20-41	30-56
Base width	58	66	76
Base length	89	100	112
Tray inside width (widest p	oint) 57	57	75
Tray inside length	44	44	44
Max. working load (kg)	45	68	91



					E930 small		E940 medium		E950 large	
			Prone Stander includes a pair of trunk		E931		E941		E951	
-0	laterals, hip str Specify colo		, padding and footboard.	CH00SE		RED 🗌		RED □		RED □
		Specify colorifica of state		m		BLU 🗌	BLU 🗆		BLU 🗆	
Optio	Additional trunk laterals (pair)				N/A		E614		E614	
nal acc	Knee laterals (pair)				E942		E942		E942	
Optional accessories	Additional wide strap				N/A		E946		E956	
Ñ	Tray (see dimension chart for size)				E947		E947		E957	
	Round abductor				E933		E933		E933	
	Collar for round abductor				E934		E934		E934	
	Round abductor with hip stabilizer				E655		E655		E655	
	djustable abduction wedge			N/A		E945		E945		
	Adjustable abduction wedge with hip stabilizer				N/A		E605		E605	
	Sandals (pair)	Ω	Small		E422		E422		N/A	
		снооѕе	Medium		E433		E433		E433	
		m	Large		N/A		N/A		E432	
	Wedges (pair, requires sandals)	요	Small		E424		E424		N/A	
		снооѕе	Medium		E435		E435		E435	
		m	Large		N/A		N/A		E434	
	Hand Anchor (For more hand/arm support options see pp. 98–101.)				K821		K821		K821	
Tota	tal			E930		E940		E950		