How to buy

- **1. Measure** your client to determine the size you need. (See dimension chart for help.)
- **2. Select** the appropriate size and color.
- **3. Choose** the accessories you need.



The **conversion kit** converts our smallest stander, the K110, to the next size up, the K130, for less than one quarter the cost of a new stander. Request from customer service.

Conversion kit K116

User dimensions (cm)	K110	K130	K150	K170
	mini	small	med.	large
Height	64-81	76-102	94-132	122-165



Key user dimension: height

Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.

Important: User's weight must not exceed the maximum working load.

Item dimensions (cm)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	64 x 64	64 x 65	81 x 74	99 x 85
Overall height	57-71	60-76	79-99	95-124
Width without big wheels	50	53	61	72
Wheel diameter	51	51	69	89
Width of body support	18-25	18-25	23-30	25-36
Height of body support	42-56	56-74	74-97	91-122
Item weight (kg)	13	14	18	29
Max. working load (kg)	23	23	54	79

